




Electrode to Connector Map

Assembly type	Any 'C' style 16-channel probe on C16 package (standard).	
Version	V6. July 2007	
Connector	Omnetics nano NPD series (A7695-001)	
Mating Connector	Omnetics nano NPD series (A7519-001)	
Compatible Headstages	TDT	RA16CH, LP16CH
	Multi-Channel Systems	Two MPA8I
	Plexon	HST/16o25-18P, HST/16V-G20
	Neuralynx	HS-16
	Alpha Omega	16 mini, TeleSpike
	DataWave Technologies	DW-T-HS16, DW-W-HS16, DW-M-HS16
Compatible Adaptors	TDT	AC-CH
Reference site	Present on some designs. Please check the catalog CAD drawings.	
Note	Electrode sites face down (as shown in the picture).	

CONNECTOR PIN OUT



ELECTRODE SITE MAP

