


Electrode to Connector Map		
Assembly type	Any 'A' style 16-channel probe on F16 package.	
Version	V1. July 2007	
Connector	Omnetics nano NPD series (A7996-001)	
Mating Connector	Omnetics nano NPD series (A7519-001)	
Compatible Headstages	TDT	RA16CH, LP16CH
	Multi-Channel Systems	Two MPA8I
	Plexon	HST/16o25-18P, HST/16V-G20
	Neuralynx	HS-16
	Alpha Omega	16 mini, TeleSpike
	DataWave Technologies	DW-T-HS16, DW-W-HS16, DW-M-HS16
Compatible Adaptors	TDT	AC-CH
Reference site	Present on limited 16-channel designs. Please check the catalog CAD drawings.	
Note	Electrode sites face up (as shown in the picture).	

CONNECTOR PIN OUT



ELECTRODE SITE MAP

