


## Electrode to Connector Map

Assembly type	Any 'A' style 16-channel probe on CM16 package.	
Version	V6. July 2007	
Connector	Omnetics nano NPD series (A7996-001)	
Mating Connector	Omnetics nano NPD series (A7519-001)	
Compatible Headstages	TDT	RA16CH, LP16CH
	Multi-Channel Systems	Two MPA8I
	Plexon	HST/16o25-18P, HST/16V-G20
	Neuralynx	HS-16
	Alpha Omega	16 mini, TeleSpike
	DataWave Technologies	DW-T-HS16, DW-W-HS16, DW-M-HS16
Compatible Adaptors	TDT	AC-CH
Reference site	N/A	
Note		

## CONNECTOR PIN OUT



## ELECTRODE SITE MAP

