

A64			
Package	SmartLink	Package	SmartLink
1	48	33	8
2	52	34	14
3	53	35	12
4	56	36	27
5	57	37	28
6	58	38	30
7	50	39	26
8	61	40	24
9	54	41	31
10	63	42	22
11	59	43	29
12	47	44	20
13	60	45	25
14	45	46	18
15	62	47	23
16	43	48	17
17	46	49	21
18	41	50	1
19	44	51	19
20	39	52	3
21	42	53	16
22	34	54	4
23	40	55	0
24	32	56	9
25	38	57	2
26	36	58	13
27	33	59	5
28	35	60	6
29	37	61	7
30	51	62	10
31	49	63	11
32	55	64	15

A32	
Package	SmartLink
21	0
17	1
23	2
18	3
22	4
19	5
25	6
20	7
27	8
24	9
29	10
26	11
28	12
31	13
30	14
32	15
11	16
9	17
7	18
5	19
1	20
3	21
4	22
2	23
13	24
6	25
14	26
8	27
15	28
10	29
16	30
12	31

A16	
Package	SmartLink
9	0
10	1
11	2
12	3
13	4
14	5
15	6
16	7
1	8
2	9
3	10
4	11
5	12
6	13
7	14
8	15