

CM16	
Package	SmartLink
1	14
2	12
3	10
4	8
5	15
6	13
7	11
8	9
9	7
10	5
11	3
12	1
13	6
14	4
15	2
16	0

H16	
Package	SmartLink
1	9
2	13
3	15
4	8
5	14
6	12
7	10
8	11
9	5
10	4
11	2
12	0
13	6
14	1
15	3
16	7

HC16	
Package	SmartLink
1	15
2	11
3	9
4	7
5	1
6	3
7	5
8	13
9	12
10	4
11	2
12	0
13	6
14	8
15	10
16	14

H16 ECoG	
Package	SmartLink
1	14
2	12
3	10
4	8
5	15
6	13
7	11
8	9
9	7
10	5
11	3
12	1
13	6
14	4
15	2
16	0