

CM32	
Package	SmartLink
1	21
2	23
3	20
4	25
5	22
6	27
7	24
8	29
9	26
10	31
11	28
12	30
13	18
14	16
15	19
16	17
17	15
18	13
19	14
20	12
21	0
22	2
23	1
24	4
25	3
26	6
27	5
28	8
29	7
30	10
31	9
32	11

H32	
Package	SmartLink
1	16
2	18
3	20
4	24
5	26
6	28
7	17
8	19
9	21
10	25
11	27
12	29
13	31
14	23
15	30
16	22
17	8
18	0
19	9
20	1
21	3
22	5
23	7
24	11
25	13
26	15
27	2
28	4
29	6
30	10
31	12
32	14

H32 ECoG	
Package	SmartLink
1	31
2	29
3	27
4	25
5	23
6	21
7	19
8	17
9	30
10	28
11	26
12	24
13	22
14	20
15	18
16	16
17	14
18	12
19	10
20	8
21	6
22	4
23	2
24	0
25	15
26	13
27	11
28	9
29	7
30	5
31	3
32	1