

H64 1-32	
Package	SmartLink
1	48
2	49
3	50
4	51
5	52
6	53
7	54
8	55
9	56
10	57
11	58
12	59
13	60
14	61
15	62
16	63
17	33
18	35
19	37
20	41
21	43
22	45
23	46
24	44
25	42
26	38
27	36
28	34
29	32
30	40
31	47
32	39

H64 33-64	
Package	SmartLink
33	24
34	16
35	23
36	31
37	29
38	27
39	25
40	21
41	19
42	17
43	18
44	20
45	22
46	26
47	28
48	30
49	0
50	1
51	2
52	3
53	4
54	5
55	6
56	7
57	8
58	9
59	10
60	11
61	12
62	13
63	14
64	15

H64LP 1-32	
Package	SmartLink
1	48
2	49
3	50
4	51
5	52
6	53
7	54
8	55
9	56
10	57
11	58
12	59
13	60
14	61
15	62
16	63
17	33
18	32
19	35
20	37
21	36
22	39
23	41
24	40
25	43
26	45
27	44
28	47
29	46
30	42
31	38
32	34

H64LP 33-64	
Package	SmartLink
33	29
34	25
35	21
36	17
37	16
38	19
39	18
40	20
41	23
42	22
43	24
44	27
45	26
46	28
47	31
48	30
49	0
50	1
51	2
52	3
53	4
54	5
55	6
56	7
57	8
58	9
59	10
60	11
61	12
62	13
63	14
64	15